

GENUS UNI

Unicompartmental knee
replacement

ADLER[®]
ORTHO



GENUS UNI Respect of the patients' anatomy.

- Main indications for unicompartmental knee arthroplasty
- Primary osteoarthritis of one compartment (medial or lateral).
- Good functionality of cruciate and collateral ligaments
- Not more than 15° of varus / valgus deformity and 10° of flexion deformity
- Absence of any degenerative disease of a progressive nature (eg Rheumatoid Arthritis)




Genus UNI has been designed according to the following principles:

- Respect for the morphology and anatomy of the patient.
- Scrupulous respect for the ligamentous apparatus.
- Minimum bone sacrifice
- Logical, reproducible and minimally invasive surgical Technique.

Instruments

The Instrument set allows the accurate reconstruction of the joint thanks to a special system of spacers.

Millimetric spacers can be added to the handle in order to check extension and flexion spaces and calibrate the bone resection depending on the cartilage wear of the patient.




Millimetric spacer to compensate cartilage wear in extension.

Handle for the control of spaces in flexion / extension

Millimetric spacer to determine the space in flexion.



Check of space in flexion



Check of space in extension


In this way we can avoid overcorrection, which is responsible for the failure of many unicompartmental knee replacements.

GENUS UNI Complete and Innovative.

Options.


Two possible primary fixation

Both tibial and femoral components are available in cemented or cementless version.



Genus UNI cementless femur.


Manufactured using the powders technology, it features an HA coated three-dimensional Co-Por® monolithic surface.



Magnification of Ti-Por® Surface.

Genus UNI cemented femur.

It has been optimized with discharging compression of the cement in the anterior and posterior distal parts.



Genus UNI – cementless tibial component.

Manufactured using the powders technology. The bone fixation is guaranteed by the HAP coated three-dimensional Ti-Por® monolithic surface.




Genus UNI – Cemented tibial component.

There are anterior discharges for cementing in order to reduce the risk of leaving residual cement in the posterior part.

An All Poly tibial component is also available.

Six femoral and tibial sizes.

Each tibial and femoral component has been designed to improve the fit with the internal or external compartment that it will replace.





Five tibial thickness for each size.

Both "Metal-Back" and "All Poly" versions are available in 5 tibial thicknesses (9mm, 10mm, 11mm, 12mm and 14mm) to achieve a millimetric optimization of the joint stability.




Femoral Component.

- Bone saving
- Balance of ligaments
- Flexion
- Anatomical
- Safety




Bone saving

6mm constant thickness of the implant.
Better control of space in flexion and extension with minimal bone resection.



Balance of ligaments.

Rotation center at constant radius of flexion from **0° to 90°**




Flexion

The **3°** closure of the posterior condyle increases the range of motion and the stability of the system.




Anatomical

External distal apex blunted. Optimizes M/L positioning without risk of protrusion.



The prosthesis follows the anatomy of the medial and lateral condyles respecting the mechanical axis of the system.



Safety

Outer edge with softened radius. Minimizes the risk of damage on polyethylene and tip effect.


Tibial Component.

- Anatomical
- Safety
- Stability
- Enhanced insert



Anatomical

The shape and dimensions of the tibial component have been designed to adapt to both the medial and lateral compartment, maximizing the contact with the peripheral cortical bone.





Safety

Dovetail implant blocking system to minimise micro-movements. Mirror finish to reduce the risk of generating polyethylene debris to a minimum.




Stability

Stabilization keel to contrast varus / valgus stress.



Anchoring pin for primary stabilization of the tibial tray. Can be used as a path for the insertion of a fixing screw (cementless option)



Enhanced Insert.

Anterior chamfer to avoid potential patellar conflicts.



Outer edge rounded to cover the Metal Back, to avoid potential conflicts between metal and soft parts.



Flat articular surface to facilitate a linear contact with the condyle and better distribute the loading forces.



Adler Ortho SpA
Via dell'Innovazione 9
20032 Cormanico, Italy
Tel +39 02 6154371
Fax +39 02 615437222
www.adlerortho.com

Adler Ortho UK
The Stables
Tarvin Road
Frodsham - Cheshire -
WA6 6XN
Tel: +44151 329 3372



B0022-00

06-2019

Femoral Components

| SIZE | Ref. Codes Cemented | | Ref.Codes Cementless (*) | |
|------|---------------------|---------|--------------------------|---------|
| | RM/LL | LM/RL | RM/LL | LM/RL |
| 1 | 1182001 | 1182011 | 1115001 | 1115011 |
| 2 | 1182002 | 1182012 | 1115002 | 1115012 |
| 3 | 1182003 | 1182013 | 1115003 | 1115013 |
| 4 | 1182004 | 1182014 | 1115004 | 1115014 |
| 5 | 1182005 | 1182015 | 1115005 | 1115015 |
| 6 | 1182006 | 1182016 | 1115006 | 1115016 |

Tibial Bases

| SIZE | Ref. Codes Cemented | | Ref.Codes Cementless (*) | |
|------|---------------------|---------|--------------------------|---------|
| | RM/LL | LM/RL | RM/LL | LM/RL |
| 1 | 1784001 | 1784011 | 1735001 | 1735011 |
| 2 | 1784002 | 1784012 | 1735002 | 1735012 |
| 3 | 1784003 | 1784013 | 1735003 | 1735013 |
| 4 | 1784004 | 1784014 | 1735004 | 1735014 |
| 5 | 1784005 | 1784015 | 1735005 | 1735015 |
| 6 | 1784006 | 1784016 | 1735006 | 1735016 |

(*) Ingrowth Surface Ti-Po+ HA Coating

Tibial Screws

| CODE | DESCRIPTION | length(mm) |
|---------|------------------------|------------|
| 1710020 | Genus Uni Tibial Screw | 20 |
| 1710025 | Genus Uni Tibial Screw | 25 |
| 1710030 | Genus Uni Tibial Screw | 30 |
| 1710035 | Genus Uni Tibial Screw | 35 |
| 1710040 | Genus Uni Tibial Screw | 40 |
| 1710045 | Genus Uni Tibial Screw | 45 |
| 1710050 | Genus Uni Tibial Screw | 50 |

(*) Available on request

All Poly tibial components

| SIZE | Ref.CodeS H9 | | Ref.CodeS H10 | | Ref.CodeS H11 | | Ref.CodeS H12 | | Ref.CodeS H14 | |
|------|--------------|---------|---------------|---------|---------------|---------|---------------|---------|---------------|---------|
| | RM/LL | LM/RL | RM/LL | LM/RL | RM/LL | LM/RL | RM/LL | LM/RL | RM/LL | LM/RL |
| 1 | 1740901 | 1740911 | 1741001 | 1741011 | 1741101 | 1741111 | 1741201 | 1741211 | 1741401 | 1741411 |
| 2 | 1740902 | 1740912 | 1741002 | 1741012 | 1741102 | 1741112 | 1741202 | 1741212 | 1741402 | 1741412 |
| 3 | 1740903 | 1740913 | 1741003 | 1741013 | 1741103 | 1741113 | 1741203 | 1741213 | 1741403 | 1741413 |
| 4 | 1740904 | 1740914 | 1741004 | 1741014 | 1741104 | 1741114 | 1741204 | 1741214 | 1741404 | 1741414 |
| 5 | 1740905 | 1740915 | 1741005 | 1741015 | 1741105 | 1741115 | 1741205 | 1741215 | 1741405 | 1741415 |
| 6 | 1740906 | 1740916 | 1741006 | 1741016 | 1741106 | 1741116 | 1741206 | 1741216 | 1741406 | 1741416 |

Inserts

| SIZE | Ref.CodeS H9 | | Ref.CodeS H10 | | Ref.CodeS H11 | | Ref.CodeS H12 | | Ref.CodeS H14 | |
|------|--------------|---------|---------------|---------|---------------|---------|---------------|---------|---------------|---------|
| | RM/LL | LM/RL | RM/LL | LM/RL | RM/LL | LM/RL | RM/LL | LM/RL | RM/LL | LM/RL |
| 1 | 1850901 | 1850911 | 1851001 | 1851011 | 1851101 | 1851111 | 1851201 | 1851211 | 1851401 | 1851411 |
| 2 | 1850902 | 1850912 | 1851002 | 1851012 | 1851102 | 1851112 | 1851202 | 1851212 | 1851402 | 1851412 |
| 3 | 1850903 | 1850913 | 1851003 | 1851013 | 1851103 | 1851113 | 1851203 | 1851213 | 1851403 | 1851413 |
| 4 | 1850904 | 1850914 | 1851004 | 1851014 | 1851104 | 1851114 | 1851204 | 1851214 | 1851404 | 1851414 |
| 5 | 1850905 | 1850915 | 1851005 | 1851015 | 1851105 | 1851115 | 1851205 | 1851215 | 1851405 | 1851415 |
| 6 | 1850906 | 1850916 | 1851006 | 1851016 | 1851106 | 1851116 | 1851206 | 1851216 | 1851406 | 1851416 |